

Western Blotting Protocol

Product Information

Product Name: Catalog Number: anti-Caspase-1 (p20) (human), mAb (Bally-1) AG-20B-0048

Caspase-1 (human) Detection by Western Blotting:

1. Priming:

THP1 cells are cultured in RPMI 1640 media, supplemented with 10% FBS. Cells are treated for 3h with phorbol-12-myristate-13-acetate (PMA) (100nM) (#AG-CN2-0010) or LPS (50ng/ml) (#AG-CU1-0001).

2. Inflammasome activation:

PMA or LPS-differentiated THP1 cells are plated in 12-well plates and treated for 6h with monosodium urate (MSU) (150µg/ml) (#AG-CR1-3951) or for 30min with Nigericin (5µM) (#AG-CN2-0020).

During this step, pro-Caspase-1 is cleaved into p10 and p20 fragments and secreted.

3. Western blotting:

Supernatants or Cell extracts from activated or non-activated cells are processed for immunoblotting, separated on a 15 % SDS-PAGE gel and electrotransferred to nitrocellulose or PVDF membrane.

 Human Caspase-1 (p20) is detected by Caspase-1 (p20) (human), mAb (Bally-1) (1µg/ml) (#AG-20B-0048)

23-10-2018/AT

WARNING: Intended for research use only. This product is not intended or approved for human, diagnostics, therapeutic or veterinary use. Use of this product for human or animal testing is extremely hazardous and may result in disease, severe injury, or death. MATERIAL SAFETY DATA: Review the complete Material Safety Data Sheet before use.

Adipogen AG Schützenstrasse 12 4410 Liestal Switzerland TEL +41-61-926-60-40 FAX +41-61-926-60-49 info@adipogen.com Adipogen Corp. 9853 Pacific Heights Blvd., Suite L San Diego CA 92121-4721 USA TEL (858) 457-8383 FAX (858) 457-8484 info-us@adipogen.com

www.adipogen.com